

## QC and cleanup

anomalous  
signal removal



Open all files in FlowJo.

Ensure compensation matrix is correct

**Exclude** irrelevant events by gating out dead cells and aggregates

**Export** live single cells with only **compensated channels**  
into new FCS files

open exported files in a new FJ workspace

### Annotate channels

#### Panel A


CD45	BUV805
CD19	missing
IgM	missing
IgD	missing
CD4	BUV661
Foxp3	PE

#### Panel B

CD45	missing
CD19	missing
IgM	missing
IgD	missing
CD4	BUV805
Foxp3	PE

#### Panel C

CD45	missing
CD19	BUV570
CD3/IgD	BUV395
CD4	BUV805
Foxp3/IgM	PE



create panel-specific  
**gating strategies** and  
annotate  
**equivalent populations**  
using one common name

{ } Blood

{ } BM

{ } Ileum

{ } Liver

{ } Lung

{ } Donor 390C

{ } Donor 403C

{ } Donor 412C

{ } Donor 423C

{ } Donor 428C

assign all samples to  
**groups** to ease data  
filtering in FA

import into FA  
explore

add new  
populations  
to FJ

# FLOWAtlas

**filter** and  
**colour** by conditions  
and groups

### Populations

- B cells
- CD4 | Naive
- CD4 | EM
- CD8 | EMRA
- CD4 | CM
- CD8


### Conditions

- BM
- Liver
- Ileum


### Groups

- 390C
- 403C
- 412C

**zoom, pan, draw,**  
**discover populations**



### differential analysis of expression



### differential analysis of population abundances

